



Your Partner In Protection™

# Surgical Gloves

For more details on specifications, please contact us at [sales@wrpworld.com](mailto:sales@wrpworld.com)

## Main Features

- Exclusive Damp-Hand Donning (DHD™) coating provides smooth and easy donning
- Custom formulation is both strong and comfortable
- "Skin Grip" outer surface is the perfect balance of tack and smoothness
- Polyisoprene film provides the gold standard in latex-free comfort
- Nitrile film provides a very economical latex-free option
- Powdered and Powder-Free products available
- Long and Short sterile pouch options
- Radiation Protection (RadiaXon only)
- Latex-free options for those sensitive to latex proteins
- Beaded Cuff

# Surgical Gloves



Your Partner In Protection™

## Surgical Powdered Latex

( Article No : WRP SG NR PP 12 )

Style No	Material	Coating	Weight per pair, size 7.5, g (min)	Sizes	Palm Thickness mm(min)	Length, mm	Colours
WRP SG NR PP 12 - 12	Natural rubber	Powder	22	5.5 - 9.0	0.17	290	White
WRP SG NR PP 12 - 9	Natural rubber	Powder	16	5.5 - 9.0	0.11	290	White

## Surgical Powder Free Latex

( Article No : WRP SG NR PF 125 )

Style No	Material	Coating	Weight per pair, size 7.5, g (min)	Sizes	Palm Thickness mm(min)	Length, mm	Colours
WRP SG NR PF 125 - PCP	Natural rubber	Polymer Coated	22	5.5 - 9.0	0.13	290	White
WRP SG NR PF 125 - PTN DHD	Natural rubber	DHD	22	5.5 - 9.0	0.13	290	White
WRP SG NR PF 125 - Micro DHD	Natural rubber	DHD	20	5.5 - 9.0	0.11	290	Brown
WRP SG NR PF 125 - EP DHD	Natural rubber	DHD	25	5.5 - 9.0	0.15	290	Brown
WRP SG NR PF 125 - SN DHD	Natural rubber	DHD	20	5.5 - 9.0	0.15	290	Green

## Surgical Powder Free Synthetic

( Article No : WRP SG PI PF 135 | WRP SG CR PF 145 | WRP SG NBR PF 75 )

Style No	Material	Coating	Weight per pair, size 7.5, g (min)	Sizes	Palm Thickness mm(min)	Length, mm	Colours
WRP SG PI PF 135 - PI DHD	Polyisoprene	DHD	22	5.5 - 9.0	0.13	290	White
WRP SG PI PF 135 - Micro DHD	Polyisoprene	DHD	20	5.5 - 9.0	0.11	290	White
WRP SG PI PF - 135 EP DHD	Polyisoprene	DHD	25	5.5 - 9.0	0.15	290	Brown
WRP SG PI PF 135 - SN DHD	Polyisoprene	DHD	20	5.5 - 9.0	0.11	290	Green
WRP SG CR PF 145 - CR DHD	Chloroprene	DHD	22	5.5 - 9.0	0.13	290	Brown
WRP SG NBR PF 75 - EPIC N	Nitrile	Nil	14	5.5 - 9.0	0.14	290	White

## Radiation Protection - Radiaxon

Article No	Material	Coating	Weight per pair, size 7.5, g (min)	Sizes	Fingertip Thickness mm(min)	Length, mm	Colours
WRP RPG NR	Natural rubber	Nil	88.0	6.0 - 9.0	0.30	275	Brown
WRP RPG PI	Polyisoprene	Nil	98.0	6.0 - 9.0	0.30	275	Black

### Quality Assurance:

US FDA Quality System Regulation (QSR), ISO9001, BS EN ISO 13485, MDSAP

COMPLIANCE & TESTED IN ACCORDANCE  
 ASTM : D 6319- 10 (2015) , EN 455 : Part 1,2,3,4 , EU CE Regulation 2016/ 425, FDA 510 K, EN 420 : 2003 + A1 2009 , EN ISO 374-1: 2016, EN 374-2: 2014, EN 374-4: 2013, EN ISO 374-5: 2016, EN 16523: 2015, EU 10/2011



Please contact us for further enquiry.

**WRP Asia Pacific Sdn Bhd** <sup>(147817-V)</sup>

Lot 1, Jalan 3, Kawasan Perusahaan Bandar Baru Salak Tinggi, 43900 Sepang, Selangor Darul Ehsan, Malaysia

Tel : +60-3-8706 1486

Fax : +60-3-8706 1557

Email : [customer\\_wrp@wrpworld.com](mailto:customer_wrp@wrpworld.com)  
[sales@wrpworld.com](mailto:sales@wrpworld.com)